

BC 4000 Benchtop centrifuge

Take your lab's performance to the peak



The fastest way to better results

BC 4000 Benchtop centrifuge



Features

1. Widely used in laboratory, hospital and etc...
2. Brushless motor for model BC4000, free maintenance, no powder pollution, quick in speed up and down.
3. Range of speed from 0 to 4000rpm, smooth in operation, low noise and small vibration.
4. Microcomputer control system, digital display the RCF, time and speed.
5. Electric cover lock, compact design, super speed and imbalance protection.
6. With over speed and imbalance protection.

BC4000 Technical Parameter

Max. Speed	4000rpm
Max. RCF	1880Xg
Max. Capacity	8X20/15/10 ml
Time Range	0~99min
RPM/RCF Convert	Yes
Noise (dB)	≤ 55
Temperature	Normal
Acc/Dcc	10 Kinds
Speed Accuracy	±10r/min
Temperature Accuracy	/
Voltage (V/Hz)	AC 220V/110V 50HZ/60HZ
Size (W x D x Hmm)	310×270×220mm
Net Weight (Kg)	6KG
Certificates	CE,ISO & Calibration report are available

Matched Rotors for BC4000

Order No.	Rotor	Max Speed (rpm)	Max Volume (ml)	Max.RCF (×g)
BC19001	Angle Rotor	4000	8×15ml	1880
BC19002	Angle Rotor	4000	6×10ml	1800
BC19003	Angle Rotor	4000	12×7/5ml	1720
BC19004	Angle Rotor	4000	8×15ml	1800
BC19005	Angle Rotor	4000	8×7/5ml	1800

Inzek B.V.

Vissenstraat 32 | 7324 AL Apeldoorn | The Netherlands | info@biozek.com | www.biozek.com